



AM Multiplexed ctDNA Standard

NAME AND INTENDED USE

The AM Multiplexed ctDNA Standard is intended for use with molecular diagnostic assays that are designed to identify somatic mutations in DNA from human samples. The assays can be either amplicon- or capture-based. The control is intended to provide a true-human plasma-based standard for validation of the assay, assessing day-to-day test variation and may help in identifying increases in random or systematic error, such as reagent lot changes, operator-based deviations, and instrument malfunction. This product is for Research Use Only.

CATALOG NUMBER

60100103, Multiplexed ctDNA 3% AF Standard (300ng in TE).

60100104, Multiplexed ctDNA 3% AF Standard (25ng/ml, 5ml, in Plasma).

SUMMARY

The AM Multiplexed ctDNA Standard contains well defined wildtype genomic DNA (gDNA) nucleosomally fragmented to around 145 bp. It carries 23 mutations formulated at around 3% allele frequency. It also contains heterozygous mutations, BRCA1 (P871L) and BRCA2 (N372H). The mixture of DNA is spiked either in DNA-depleted human plasma or in TE-based buffer. Please refer to the product Certificate of Analysis for the exact DNA concentration.

WARNINGS AND PRECAUTIONS

For Research Use Only. Not for use in diagnostic procedures. CAUTION: Handle the plasma sample as though it is capable of transmitting infectious agents.

SAFETY PRECAUTIONS

Use Centers for Disease Control and Prevention (CDC) recommended universal precautions for handling reference materials and human specimens. Do not pipette by mouth. Do not smoke, eat, or drink in areas where specimens are being handled. Clean any spillage by immediately wiping with 0.5% sodium hypochlorite solution. Dispose of all specimens and materials used in testing as though they contain infectious agents.

STORAGE INSTRUCTIONS

Store AM Multiplexed ctDNA Standard can be stored for up to 30 days at 2-8°C. For longer term storage, please keep it frozen at -20 °C or colder. Shelf life when stored under these conditions is three years from date of manufacture.

INDICATIONS OF REAGENT INSTABILITY OR DETERIORATION

AM Multiplexed ctDNA Standard is a mixture of human genomic DNA. In plasma matrix, it should appear as a yellowish or brownish liquid. In TE buffer, it should appear as a clear liquid. Alterations in this appearance may indicate instability or deterioration of the product and vials should be discarded.

HANDLING PRECAUTIONS

Do not use the product beyond the expiration date. Avoid contamination of the product when opening and closing the vial.

INSTRUCTION FOR USE

For plasma samples, after thawing, 0.5 to 5 ml of the plasma-based AM Multiplexed ctDNA Standard is extracted with a cfDNA extraction method before library preparation or direct assay. For buffer samples, the DNA may be diluted to an appropriate concentration for library preparation or direct quantitation, refer to your assay procedures in order to determine the amount of material to use.

EXPECTED RESULTS & INTERPRETATION OF RESULTS

In-house quantitation results were provided in the Certificate of Analysis (CoA). Each laboratory must establish an assay-specific expected value. When results for the product significantly deviate from values on the CoA, it may indicate unsatisfactory test performance. Possible sources of error include: deterioration of test kit reagents, operator error, faulty performance of equipment, contamination of reagents, etc.

MUTATIONS INCLUDED

| Variant Name | Cosmic ID |
|----------------------|--------------|
| KRAS G12C | COSV55497469 |
| KRAS G12D | COSV55497369 |
| KRAS G12V | COSV55497419 |
| TP53 R248Q | COSV52661580 |
| PIK3CA E545K | COSV55873239 |
| PIK3CA H1047R | COSV55873195 |
| EGFR T790M | COSV51765492 |
| EGFR C797S | COSV51766493 |
| EGFR D770_N771insG | COSV51769298 |
| EGFR L858R | COSV51765161 |
| BRAF V600E | COSV56056643 |
| NRAS Q61R | COSV54736340 |
| PTEN R130G | COSV64288384 |
| TP53 R175H | COSV52661038 |
| IDH1 R132H | COSV61615239 |
| CTNNB1 T41A | COSV62687862 |
| CTNNB1 S33C | COSV62687839 |
| FGFR3 S249C | COSV53390026 |
| APC Q1429* | COSV57322632 |
| APC R1450* | COSV57321313 |
| EGFR E746-A750del v1 | COSV51765119 |
| EGFR E746-A750del v2 | COSV51765066 |
| KIT D816V | COSV55386424 |

1. Note BRCA1 (P871L) and BRCA2 (N372H) are present as heterozygous (50% AF) mutations.