



Certificate of Analysis

Material Name:	10% BCR-ABL1 RNA
Part Number:	RF164-028, RF164-044, and RF164-045
Fill Volume:	25 ul / vial
Lot Number:	See table below
Composition and Usage:	Purified RNA containing mRNA for BCR-ABL1 p210e14a2, p210e13a2, p190e1a2 Fusion and wildtype ABL1. The product can be used directly in most PCR-based or NGS applications.
Storage:	-70°C or lower.
Date of Manufacturing:	5/12/2026
Expiration Date:	5/31/2027

Test Result (ddPCR)


Part Number	Lot Number	Analyte	ABL1 (copy/ul)	Fusion/ABL1 %	MR	Criteria (MR)	PASS/FAIL
RF164-028	604301	p210e14a2	1560	15.4%	0.81	1 +/- 0.3	PASS
RF164-044	605042	p210e13a2	1580	13.1%	0.88	1 +/- 0.3	PASS
RF164-045	604293	p190e1a2	1380	8.2%	1.09	1 +/- 0.3	PASS

Note: The "Fusion/ABL1 %" value for p210 e14a2 is traceable to the WHO International Scale (IS) via calibration to the NIBSC reference material 09/138.

Approval

Name: Yabin Lu

Signature:

DocuSigned by:

A68122C2A4CC4F1...

Date: 5/13/2026

The material is for research use only.